




Microbiology

TOTAL URINE



**REAL URINE WITH BIOCHEMISTRY AND CYTOBACTERIOLOGICAL
PARAMETERS PERSONALIZED TO YOUR PREFERENCES**

Patented technique using real urine

Produced in France with ISO 13485 certification

6-month-long shelf life

TOTAL URINE

The product is composed of two vials to be mixed:


- ◆ A freeze dried vial with bacteriological parameters
- ◆ A liquid vial with real urine containing cytological and biochemistry parameters

▶ FREEZE DRIED VIAL WITH BACTERIOLOGICAL PARAMETERS

- ◇ **Crystals:** Oxalate, Struvite, Cholesterol
- ◇ **Bacteria:** Strains, Enumeration, Phenotypic resistance


▶ LIQUID VIAL WITH CYTOLOGICAL AND BIOCHEMISTRY PARAMETERS

- ◇ **Red Blood Cells and/or White Blood Cells:** 0 to 10^6 mL
- ◇ Select biochemistry parameters and their concentration according to the following list:
**Amylase - Calcium - Chlorine - Creatinine - Glucose - Magnesium - Microalbumin
Phosphate - Potassium - Protein - Sodium - Urea - Uric acid**



Struvite and Cholesterol crystals

OUR EQA SAMPLES ARE FULLY CUSTOMISABLE IN TERMS OF:



CONTACT INFORMATION

ELITech Microbio
Parc d'Activités, Allée d'Athènes
83870 Signes – France
Tel: +33 (0)4 94 88 55 00
Fax: +33 (0)4 94 32 82 61

www.elitechgroup.com
em-marketing@elitechgroup.com



: ELITech Microbio (France)

ELITechGroup
A BRUKER COMPANY